

Version 1a Last updated 6 December 2018

ab233471 Lipid Peroxidation (MDA) Assay Kit (Colorimetric)

For the measurement of lipid peroxidation.

This product is for research use only and is not intended for diagnostic use.

Table of Contents

1. Overview	1
2. Materials Supplied and Storage	2
3. Materials Required, Not Supplied	3
4. General guidelines, precautions, and troubleshooting	4
5. Reagent Preparation	5
6. Standard Preparation	6
7. Sample Preparation	6
8. Assay Procedure	7
9. Data Analysis	9
10. FAQs / Troubleshooting	10
11. Typical Data	11
12. Notes	12

1. Overview

Lipid Peroxidation (MDA) Assay Kit (Colorimetric) (ab233471) offers the most rapid and convenient method to measure malondialdehyde (MDA) without TBARS heating steps. MDA Color Reagent reacts with MDA to generate a blue color product which is measured at 695 nm with absorbance microplate readers. This assay is very fast and specific to MDA with little interference from other aldehydes.

Prepare serially diluted MDA standards and test samples.



Add 10 μ l of MDA Color Reagent stock solution into each well of MDA standard.



Incubate at room temperature for 10-30 minutes.



Add 40 μ l of Reaction Solution and incubate at room temperature for 30-60 minutes.



Monitor absorbance increase at 695 nm.

2. Materials Supplied and Storage

Store kit at -20°C immediately upon receipt. Kit can be stored for 1 year from receipt, if components have not been reconstituted.

Aliquot components in working volumes before storing at the recommended temperature.

Avoid repeated freeze-thaws of reagents.

Item	Quantity	Storage temperature (before prep)	Storage temperature (after prep)
MDA Color Reagent	1 vial	-20°C	-20°C
Dilution Buffer	10 mL	-20°C	-20°C
MDA Standard	1 vial	-20°C	-20°C
Reaction Solution	10 mL	-20°C	-20°C

3. Materials Required, Not Supplied

These materials are not included in the kit, but will be required to successfully perform this assay:

- Ultra pure and sterile water (ddH₂O).
- 96-well clear bottom microplates.
- Microplate reader capable of reading absorbance at 695 nm.

4. General guidelines, precautions, and troubleshooting

Please observe safe laboratory practice and consult the safety datasheet.

For general guidelines, precautions, limitations on the use of our assay kits and general assay troubleshooting tips, particularly for first time users, please consult our guide:

www.abcam.com/assaykitguidelines

For typical data produced using the assay, please see the assay kit datasheet on our website.

5. Reagent Preparation

Briefly centrifuge small vials at low speed prior to opening.

5.1 MDA Color Reagent

Ready to use as supplied.

ΔNote: Aliquot into single use size and store at -20°C protected from light.

5.2 Dilution Buffer

Ready to use as supplied.

5.3 MDA Standard

Add 100 μ L of ddH₂O into MDA Standard vial to make 100 mM MDA stock solution.

5.4 Reaction Solution

Ready to use as supplied.

6. Standard Preparation

Prepare serially diluted MDA standards

1. Prepare 100 mM MDA stock solution as described above (Reagent Preparation 5.3).
2. Prepare MDA standard dilutions: add 4 μL of 100 mM MDA stock solution (from step 6.1.1.) into 996 μL of Dilution Buffer to obtain 400 μM MDA solution.
3. Perform 1:2 serial dilutions in Dilution Buffer to obtain 200, 100, 50, 25, 12.5, 6.25 and 0 μM serially diluted MDA standards.

7. Sample Preparation

1. Prepare 1:2 serially diluted Test Samples in Dilution Buffer.

8. Assay Procedure

- Equilibrate all materials and prepared reagents to room temperature just prior to use and gently agitate.
- Assay all standards, controls and samples in duplicate.

8.1 Prepare assay plate

Add 50 μl of test samples and serially diluted MDA standards into a 96-well clear bottom microplate as described in the table, below.

Blank Control	Blank Control	Test Sample	Test Sample								
MDA 1	MDA 1								
MDA 2	MDA 2	...									
MDA 3	MDA 3										
MDA 4	MDA 4										
MDA 5	MDA 5										
MDA 6	MDA 6										
MDA 7	MDA 7										

Blank Control = Dilution Buffer.

MDA 1...MDA 7 = Serially diluted MDA Standard (6.25 to 400 μM , respectively).

8.2 Run MDA assay

1. Add 10 μL MDA Color Reagent solution into each well of MDA Standard, Blank Control, and Test Samples.

Δ Note: For a 384-well plate, add 25 μL of MDA Standard /Test Sample and 5 μL of MDA Color Reagent stock solution into each well.

2. Incubate the reaction mixture at room temperature for 10-30 minutes.
3. Add 40 μL of Reaction Solution to make the total assay volume to 100 μL /well.

ΔNote: For a 384-well plate, add 20 μL of Reaction Solution to make the total assay volume 50 μL /well.

4. Incubate the final reaction mixture at room temperature for 30-60 minutes.
5. Monitor absorbance increase with an absorbance plate reader with path-check correction at OD of 695~700 nm.

9. Data Analysis

The absorbance reading in blank wells (with Dilution Buffer only) is used as a control, and is subtracted from the value of those wells with the standards and test samples.

Calculate Test Sample MDA concentrations by reference to the MDA Standard Curve correcting for the Test Sample dilution.

10.FAQs / Troubleshooting

General troubleshooting points are found at www.abcam.com/assaykitguidelines.

11. Typical Data

Data provided for demonstration purposes only.




Figure 1. MDA dose response was measured using Lipid Peroxidation (MDA) Assay Kit (Colorimetric)(ab233471) in a 96-well clear bottom microplate. MDA standard curve (left). Signal comparison of MDA, Formaldehyde and Glyceraldehyde (right).

12. Notes

Technical Support

Copyright © 2012-2017 Abcam. All Rights Reserved. The Abcam logo is a registered trademark. All information / detail is correct at time of going to print.

Austria

wissenschaftlicherdienst@abcam.com | 019-288-259

France

supportscientifique@abcam.com | 01.46.94.62.96

Germany

wissenschaftlicherdienst@abcam.com | 030-896-779-154

Spain

soportecientifico@abcam.com | 91-114-65-60

Switzerland

technical@abcam.com

Deutsch: 043-501-64-24 | Français: 061-500-05-30

UK, EU and ROW

technical@abcam.com | +44(0)1223-696000

Canada

ca.technical@abcam.com | 877-749-8807

US and Latin America

us.technical@abcam.com | 888-772-2226

Asia Pacific

hk.technical@abcam.com | (852) 2603-6823

China

cn.technical@abcam.com | 400 921 0189 / +86 21 2070 0500

Japan

technical@abcam.co.jp | +81-(0)3-6231-0940

Singapore

sg.technical@abcam.com | 800 188-5244

Australia

au.technical@abcam.com | +61-(0)3-8652-1450

New Zealand

nz.technical@abc.com | +64-(0)9-909-7829