

Version 1 Last updated 31 October 2018

ab241011 Alginate Hydrogel Kit for 3D Cell Culture

For spheroid formation assays in adherent and suspension cells.

This product is for research use only and is not intended for diagnostic use.

Table of Contents

1. Overview	3
2. Protocol Summary	3
3. General guidelines, precautions, and troubleshooting	4
4. Materials Supplied, and Storage and Stability	6
5. Materials Required, Not Supplied	6
6. Reagent Preparation	7
7. Assay Procedure	8
8. Typical Data	10
9. Notes	11

1. Overview

Alginate Hydrogel Kit for 3D Cell Culture (ab241011) offers 3D cell culture matrices to meet the needs and requirements of various research fields. This kits provide a standardized, yet user friendly and adaptable to high-throughput strategy for setting up spheroid formations, 3D cell cultures and pharmacological studies.

2. Protocol Summary

Prepare reagents and grow cells in appropriate culture conditions.



Harvest, centrifuge and count the number of cells. Resuspend cells in media at the concentration of 2×10^6 cells/ml.



Add 500 μ L of cells to 4.5 mL of Alginate Hydrogel Matrix at RT. Mix and add 50 μ L of cell mixture.



To solidify the matrix, add 250 μ L of ice-cold Cross-linking Solution to each well and incubate at RT for 5-10 min.



After solidification, remove all liquid and Alginate Hydrogel remains in the well. Wash three times with 200 μ L of Wash Buffer.



Add 200-250 μ L of media and allow cells to grow and form spheroids in 37°C incubator. Change media every 2-3 days.

3. General guidelines, precautions, and troubleshooting

- Please observe safe laboratory practice and consult the safety datasheet.
- For general guidelines, precautions, limitations on the use of our assay kits and general assay troubleshooting tips, particularly for first time users, please consult our guide:
www.abcam.com/assaykitguidelines
- For typical data produced using the assay, please see the assay kit datasheet on our website.

4. Materials Supplied, and Storage and Stability

- Store kit at -20°C in the dark immediately upon receipt and check below in Section 6 for storage for individual components. Kit can be stored for 1 year from receipt, if components have not been reconstituted.
- Aliquot components in working volumes before storing at the recommended temperature.

Item	Quantity	Storage condition
Alginate Hydrogel Matrix	5 mL	-20°C
Cross-linking Solution	25 mL	-20°C
Wash Buffer	100 mL	-20°C

5. Materials Required, Not Supplied

These materials are not included in the kit, but will be required to successfully perform this assay:

- Cell culture media
- 96-well plate (sterile, clear-bottom)
- Microscope
- Matrix Dissociation Buffer

6. Reagent Preparation

- Before using the kit, spin tubes and bring down all components to the bottom of tubes.
- Prepare only as much reagent as is needed on the day of the experiment.

6.1 Alginate Hydrogel Matrix:

Ready to use as supplied. Aliquot and store at -20°C . Avoid multiple freeze/thaw. Use within two months. Thaw and keep on ice before use.

6.2 Cross-linking Solution:

Ready to use as supplied. Store at -20°C . Thaw and keep at 4°C before use. Stable for six months after the first thaw.

6.3 Wash Buffer:

Ready to use as supplied. Store at -20°C or 4°C . Stable for six months after the first thaw. Bring to room temperature (RT) before use.

7. Assay Procedure

- Equilibrate all materials and prepared reagents to room temperature just prior to use and gently agitate.

7.1 Cells:

- Grow cells in appropriate media and culture conditions. Adherent cells should be cultured to ~80% confluency.
- Harvest cells and centrifuge at 1,000 x *g*, for 5 min. Resuspend the cell pellet in Wash Buffer and count the number of cells using a hemocytometer or an automated cell counter.
- Resuspend cells in 500 μL of media at the concentration of 2×10^6 cells/ml.
- For 96 well-plate, add 500 μL of resuspended cells to 4.5 mL of thawed Alginate Hydrogel Matrix at RT. Mix gently by pipetting, and add 50 μL of cell mixture to each well to get 10,000 cells per well.

Δ Note: For a scaled-down experiment to 10 wells, add 55 μL of cells in media (2×10^6 cells/ml) to 495 μL of Alginate Hydrogel Matrix. Next, add 50 μL of cell mixture to each well to get 10,000 cells per well.

7.2 Matrix Preparation:

- To solidify the matrix, add 250 μL of ice-cold Cross-linking Solution to each well. Incubate at RT for 5-10 min until matrix is formed (clear gel in liquid solution).
- After the matrix solidification, remove all liquid by pipetting, and Alginate Hydrogel remains in the well. Wash three times with 200 μL of Wash Buffer.
- Add 200-250 μL of appropriate media and allow cells to grow and form spheroids in 37°C incubator for a set amount of days depending on experimental set up. Change media every 2-3 days.

Δ Note: Cells typically form spheroids in matrix after 1 week. Matrix remains stable for up to 3 weeks in culture. Vacuum removal of buffer or media could aspirate some or the entire matrix and can cause loss of samples. Removal of media (i.e. by careful pipetting) is strongly recommended.

7.3 Matrix Dissociation (optional):

- Matrix Dissociation Buffer is not provided. Add 200-250 μL of Matrix Dissociation Solution. Incubate at RT for 5-10 minutes and then pipet up and down with 1 mL tip until matrix is dissolved. Move the cells and solution to 1.5 mL Eppendorf tubes. To neutralize the Matrix Dissociation Solution, add 1 mL of Wash Buffer to each tube and centrifuge at $1,000 \times g$, for 5 minutes. Resuspend cells in media for use in assay of interest.

8. Typical Data

Typical data provided for demonstration purposes only.




Figure 1. HepG2 and MCF7 cells in No Matrix (A) and 3D Alginate Hydrogel Matrix (B). Cells were cultured in Alginate Hydrogel Matrix for 21 days, and successfully formed spheroids. Media was changed every 2-3 days as per protocol. The Calcein AM staining (C) indicates that cell viability is not affected while culturing in matrix for a long period of time. *Note: Calcein AM is not included in the kit.*

9. Notes

Technical Support

Copyright © 2018 Abcam. All Rights Reserved. The Abcam logo is a registered trademark. All information / detail is correct at time of going to print.

Austria

wissenschaftlicherdienst@abcam.com | 019-288-259

France

supportscientifique@abcam.com | 01.46.94.62.96

Germany

wissenschaftlicherdienst@abcam.com | 030-896-779-154

Spain

soportecientifico@abcam.com | 91-114-65-60

Switzerland

technical@abcam.com

Deutsch: 043-501-64-24 | Français: 061-500-05-30

UK, EU and ROW

technical@abcam.com | +44(0)1223-696000

Canada

ca.technical@abcam.com | 877-749-8807

US and Latin America

us.technical@abcam.com | 888-772-2226

Asia Pacific

hk.technical@abcam.com | (852) 2603-6823

China

cn.technical@abcam.com | +86-21-5110-5938 | 400-628-6880

Japan

technical@abcam.co.jp | +81-(0)3-6231-0940

Singapore

sg.technical@abcam.com | 800 188-5244

Australia

au.technical@abcam.com | +61-(0)3-8652-1450

New Zealand

nz.technical@abc.com | +64-(0)9-909-7829