

Product datasheet

Fluoroshield Mounting Medium with Propidium Iodide (100ml) [ab104130](#)

Overview

Product name	Fluoroshield Mounting Medium with Propidium Iodide (100ml)
General notes	<p>Fluoroshield Mounting Medium with PI is an aqueous mounting medium for preserving fluorescence of tissue and cell smears. This unique formula prevents rapid photobleaching of FITC, Texas Red, AMCA, Cy2, Cy3, Cy5, Green fluorescent protein (GFP), tetramethyl rhodamine, Redox, Phycoerythrin (RP-E), Phyocyanin (PC), and Allophycocyanin (APC). The fluorescence is retained during prolonged storage at 4°C in the dark. This medium does not contain phenylenediamine, which destroys immunofluorescence of Cy dyes, RP-E, PC and APC. This mounting medium is fortified with PI which is a counter-stain for DNA. This product is to be used in <i>in situ</i> hybridization techniques or other methods where fluorescence of DNA staining is required. PI excites at 535nm and emits at 615nm, producing a red fluorescence. RNA is also stained with PI.</p> <p>Procedure</p> <p>Bring the vial to room temperature.</p> <p>Rinse slide to be mounted with distilled or deionized, touch the edges of slide on a paper towel to remove excess water. Place slides on a flat surface away from light.</p> <p>Turn the vial upside down and remove any air bubbles.</p> <p>Apply 3-4 drops directly of mounting medium on top of the specimen and spread out evenly by tilting the slide back and forth or spread evenly with a 0.2 ml plastic pipette tip, making sure the tissue is not touched. Excess medium can be removed by touching the edges of slide on paper towel.</p> <p>Let stand at room temperature for about 5 minutes.</p> <p>Apply cover slip, carefully avoiding air bubbles.</p> <p>The specimen is ready for visualization under microscope.</p> <p>One can seal the edges of cover slip with nail polish or any organic medium or Organo mounting medium (ab104141). If coverslip is not sealed air bubbles will appear in few days.</p> <p>For long term storage it is recommended that the slide be stored in the dark at 2-8°C. Removal of Coverslip Coverslip can be removed before sealing the edges. Soak slide in warm (37°C) water for a few minutes. Carefully and slowly move the coverslip. Soak in water an additional few minutes and remove coverslip. Rinse slide several times with warm water to remove all mounting medium. The slide can be remounted again.</p>

Properties

Storage instructions	Store at +4°C. Do Not Freeze. Store In the Dark.
Storage buffer	Preservative: 0.05% Sodium azide

Constituents: 0.0002% Propidium Iodide, 84% Water, 15% Proprietary component

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors