

Product datasheet

Anti-Pseudomonas aeruginosa serotype 5c antibody [SD6930] ab69232

★★★★★ [1 Abreviews](#) [3 References](#)

Overview

Product name	Anti-Pseudomonas aeruginosa serotype 5c antibody [SD6930]
Description	Mouse monoclonal [SD6930] to Pseudomonas aeruginosa serotype 5c
Host species	Mouse
Tested applications	Suitable for: ELISA, ICC/IF
Species reactivity	Reacts with: Pseudomonas aeruginosa
Immunogen	Tissue, cells or virus corresponding to Pseudomonas aeruginosa serotype 5c. Pseudomonas aeruginosa serotype 5c (whole organism)
General notes	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Preservative: 0.065% Sodium azide Constituent: PBS
Purity	Protein A purified
Clonality	Monoclonal
Clone number	SD6930
Myeloma	Sp2/0
Isotype	IgG1

Applications

The Abpromise guarantee

Our **Abpromise guarantee** covers the use of ab69232 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		Use at an assay dependent dilution.
ICC/IF	★★★★★ (1)	Use at an assay dependent dilution.

Target

Relevance

Pseudomonas aeruginosa is Gram-negative, aerobic, rod-shaped bacteria with unipolar motility. An opportunistic human pathogen, *P. aeruginosa* is also an opportunistic pathogen of plants. *P. aeruginosa* bacteria are clinically important because they are resistant to most antibiotics and they are capable of surviving in conditions that few other organisms can tolerate. *Pseudomonas* is often encountered in hospital and clinical work because it is a major cause of hospital acquired (nosocomal) infections. Its main targets are immunocompromised individuals, burn victims, and individuals on respirators or with indwelling catheters. Additionally, these pathogens colonize the lungs of cystic fibrosis patients. *P. aeruginosa* is often identified by its pearlescent appearance and grape-like odor in vitro. Definitive clinical identification of *P. aeruginosa* includes identifying the production of both pyocyanin and fluorescein as well as its ability to grow at 42°C. *P. aeruginosa* is capable of growth in diesel and jet fuel, where it is known as hydrocarbon utilizing microorganisms (or "HUM bugs"), causing microbial corrosion. It creates dark gellish mats sometimes improperly called "algae".

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors