abcam

Product datasheet

Recombinant Human SEC22B protein (denatured) ab134537

1 Image

Description

Product name Recombinant Human SEC22B protein (denatured)

Purity > 90 % SDS-PAGE.

Expression system Escherichia coli

Accession 075396

Protein length Protein fragment

Animal free No

Nature Recombinant

Species Human

Sequence MGSSHHHHHHSSGLVPRGSHMGSMLPLAASMQEDEQS

GRDLQQYQSQAKQ

LFRKLNEQSPTRCTLEAGAMTFHYIEQGVCYLVLCEAAFP

KKLAFAYLE

DLHSEFDEQHGKKVPTVSRPYSFIEFDTFIQKTKKLYIDSR

ARRNLGSIN

TELQDVQRIMVANIEEVLQRGEALSALDSKANNLSSLSKK

YRQDAKYLNM RSTYA

Predicted molecular weight 23 kDa including tags

Amino acids 14 to 194

Tags His tag N-Terminus

Specifications

Our $\underline{\mbox{\bf Abpromise guarantee}}$ covers the use of $\mbox{\bf ab134537}$ in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

Form Liquid

Additional notes This product was previously labelled as SEC22L1

1

Preparation and Storage

Stability and Storage

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -

80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 12.01% Urea, 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.88% Sodium

chloride

General Info

Function SNARE involved in targeting and fusion of ER-derived transport vesicles with the Golgi complex

as well as Golgi-derived retrograde transport vesicles with the ER.

Sequence similarities Belongs to the synaptobrevin family.


Contains 1 longin domain.

Contains 1 v-SNARE coiled-coil homology domain.

Cellular localization

Endoplasmic reticulum membrane. Endoplasmic reticulum-Golgi intermediate compartment membrane. Golgi apparatus, cis-Golgi network membrane. Golgi apparatus, trans-Golgi network membrane. Melanosome. Concentrated most in the intermediate compartment/cis-Golgi network and the cis-Golgi cisternae 1 and 2. Greatly reduced in concentration at the trans end of the Golgi apparatus. Identified by mass spectrometry in melanosome fractions from stage I to stage IV.

Images



15% SDS-PAGE analysis of 3 µg ab134537.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit https://www.abcam.com/abpromise or contact our technical team.

Terms and conditions

• Guarantee only valid for products bought direct from Abcam or one of our authorized distributors