

Product datasheet

Optiblot Gel Drying Solution ab119201

Overview

Product name	Optiblot Gel Drying Solution
General notes	Optiblot Gel Drying Solution is optimal for the drying of Optiblot Precast Gels for long term storage.

Optiblot gels must be used with Optiblot buffer.

Buffer recipes are available in the protocol tab.

Properties

Storage instructions	Store at +4°C.
Storage buffer	Constituent: 10% Glycerol

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors