

Product datasheet

Aqueous Mounting Medium ab128982

★★★★★ 2 Abreviews 5 References

Overview

Product name	Aqueous Mounting Medium
Tested applications	Suitable for: ICC/IF, IHC-Fr, IHC-P
General notes	<p>Aqueous mounting medium for use in fluorescence microscopy. Enhances FITC intensity and reduces fluorescence fading. Ready-to-use.</p> <p>Reagent: 6 ml</p> <p>Other products for IHC</p> <p>Find more kits and reagents for antigen retrieval, blocking, signal amplification, visualization, counterstaining, and mounting in the IHC kits and reagents guide.</p> <p>Procedure:</p> <p>FOR USE IN FLUORESCENCE MICROSCOPY</p> <p>After completion of immunostaining, remove slide from water, apply sufficient drops of Aqueous Mounting Medium upon section and mount coverglass. Coverslip should be immobilized within 30 minutes at room temperature. To avoid air bubbles within mounting medium, DO NOT SHAKE VIAL. If air bubbles are evident during coverslipping, carefully trim the vial tip with a sharp blade to create a larger opening and improve media flow. These directions are recommendations; however, they may be altered according to procedures established in-house.</p> <p>IHC Protocol Summary</p>

Properties

Form	Liquid
Storage instructions	Store at +4°C. Do not freeze.
Storage buffer	Preservative: 0.1% Sodium azide

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab128982 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC/IF		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.
IHC-P		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors