

Product datasheet

DTT Reducer (10X) - RunBlue™

ab119199

1 References

Overview

Product name	DTT Reducer (10X) - RunBlue™
Description	DTT Reducer (10X) - RunBlue™
General notes	<p>This product is manufactured by Expedeon, an Abcam company, and was previously called RunBlue™ DTT Reducer - 1ml (10X Concentrate). NXA32001 is the same as the 1 ml size.</p> <p>DTT Reducer (10X) - RunBlue™ should be added during sample preparation if you wish to perform reduced SDS-PAGE and reduce the disulfide bonds of proteins.</p> <p>RunBlue™ gels must be used with RunBlue™ buffers.</p> <p>Buffer recipes are available in the protocol tab.</p>

Properties

Form	Liquid
Storage instructions	Store at -20°C.
Storage buffer	Constituent: 31% DTT

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors