

Product datasheet

LDS Sample Buffer (4X) - RunBlue™ (Bis-Tris) ab270226

Overview

| | |
|----------------------|---|
| Product name | LDS Sample Buffer (4X) - RunBlue™ (Bis-Tris) |
| General notes | <p>LDS Sample Buffer (4X) - RunBlue™ (Bis-Tris) (ab270226) has been specifically formulated for use with Bis-Tris RunBlue™ gels to provide the best running conditions and results. The ions in the sample buffer match the gel buffer and the buffer has a higher density, making it compatible with the density of the running buffer. The buffer contains two tracking dyes. The red dye is designed to monitor the pH of the buffer. A color change from red to yellow indicates a shift in pH indicating that the buffer was poorly made or recycled too many times.</p> <p>LDS Sample Buffer (4X) - RunBlue™ (Bis-Tris) (ab270226) should be used to prepare all samples for running on Bis-Tris Precast Gels - RunBlue™.</p> <p>This product is manufactured by Expedeon, an Abcam company. It was previously called RunBlue LDS Sample Buffer - 10ml (4X Concentrate) for Bis-Tris gels. Expedeon product code NXB32010 is the same as the 10 mL size of this product.</p> <p>Please note that for TEO-Tricine Precast Gels – RunBlue™ we recommend using LDS Sample Buffer (4X) - RunBlue™ (TEO-Tricine) (ab119196).</p> |

Properties

| | |
|-----------------------------|---|
| Form | Liquid |
| Storage instructions | Shipped at Room Temperature. Store at +4°C. |
| Storage buffer | pH: 8.6 Constituents: 40% Glycerol (glycerin, glycerine), 13% Lithium dodecyl sulphate |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors