

Product datasheet

LoVo nuclear extract lysate ab14859

Overview

Product name	LoVo nuclear extract lysate
General notes	Extracts have been quality control tested by Western blot and the Electrophoretic Mobility Shift Assay (EMSA).

Properties

Mycoplasma free	Yes
Form	Liquid
Storage instructions	Aliquot and store at -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	pH: 7.90 Constituents: 0.75% Potassium chloride, 0.009% PMSF, 0.01% Magnesium chloride, 0.008% DTT, 0.48% HEPES, 20% Glycerol
Background	These epithelial like cells are used as a model for colorectal tumours as they express tumor markers and produce CEA (carcinoembryonic antigen). They are also used to investigate cytotoxicity, tumorigenicity, drug resistance, pharmacodynamics, chemotherapeutic agents as well as for antitumour testing.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors