

## Product datasheet

# Low Fluorescence Western Membrane (PVDF) ab133411

★★★★★ 1 Abreviews 4 References

### Overview

<b>Product name</b>	Low Fluorescence Western Membrane (PVDF)
<b>Tested applications</b>	<b>Suitable for:</b> WB
<b>General notes</b>	Pre-cut membranes (10 sheets, 7 x 9 cm), compatible with gels 10 x 10 cm and smaller. Note: PVDF membranes must be activated prior to use. Membrane activation protocol: 1. Pre-wet the membrane in methanol for a few moments until the membrane is completely wet. 2. Transfer the wet membrane from methanol to purified water and incubate on a rocker or shaker for at least 5 minutes. Make sure that the membrane is fully immersed in water and does not float on the water's surface. 3. Transfer the membrane from water to transfer buffer, and incubate for 5 minutes.

### Properties

<b>Form</b>	Solid
<b>Storage instructions</b>	Store at room temperature.

### Applications

**The Abpromise guarantee** Our [Abpromise guarantee](#) covers the use of ab133411 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use at an assay dependent concentration.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors