abcam

Product datasheet

PE Anti-Vimentin antibody [EPR3776] - Cytoskeleton Marker ab209446





2 References 4 Images

Overview

Product name PE Anti-Vimentin antibody [EPR3776] - Cytoskeleton Marker

Description PE Rabbit monoclonal [EPR3776] to Vimentin - Cytoskeleton Marker

Host species Rabbit

Conjugation PE. Ex: 488nm, Em: 575nm

Tested applications Suitable for: ICC/IF, Flow Cyt (Intra)

Species reactivity Reacts with: Human

Predicted to work with: Mouse, Rat, Rhesus monkey ^

Synthetic peptide. This information is proprietary to Abcam and/or its suppliers. **Immunogen**

Positive control Flow Cyt (intra): HeLa cells, wildtype HAP-1 cells ICC/IF: HeLa cells

Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit **General notes**

monoclonal antibodies. For details on our patents, please refer to **RabMAb**® **patents**.

Properties

Form Liquid

Storage instructions Shipped at 4°C. Upon delivery aliquot. Store at +4°C. Do Not Freeze. Store In the Dark.

Storage buffer pH: 7.4

> Preservative: 0.02% Sodium azide Constituents: 1% BSA, PBS

Purity Protein A purified

Clonality Monoclonal Clone number **EPR3776**

Isotype ΙgG

Applications


The Abpromise guarantee

Our <u>Abpromise guarantee</u> covers the use of ab209446 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC/IF		1/500. This product gave a positive signal in HeLa cells fixed with 4% formaldehyde (10 min) and 100% methanol (5 min)
Flow Cyt (Intra)		1/5000. ab209478 - Rabbit monoclonal lgG (Phycoerythrin), is suitable for use as an isotype control with this antibody.

Target		
Function	Vimentins are class-Ill intermediate filaments found in various non-epithelial cells, especially mesenchymal cells. Vimentin is attached to the nucleus, endoplasmic reticulum, and mitochondria, either laterally or terminally. Involved with LARP6 in the stabilization of type I collagen mRNAs for CO1A1 and CO1A2.	
Tissue specificity	Highly expressed in fibroblasts, some expression in T- and B-lymphocytes, and little or no expression in Burkitt's lymphoma cell lines. Expressed in many hormone-independent mammary carcinoma cell lines.	
Involvement in disease	Cataract 30	
Sequence similarities	Belongs to the intermediate filament family.	
Domain	The central alpha-helical coiled-coil rod region mediates elementary homodimerization. The [IL]-x-C-x-x-[DE] motif is a proposed target motif for cysteine S-nitrosylation mediated by the iNOS-S100A8/A9 transnitrosylase complex.	
Post-translational modifications	Filament disassembly during mitosis is promoted by phosphorylation at Ser-55 as well as by nestin (By similarity). One of the most prominent phosphoproteins in various cells of mesenchymorigin. Phosphorylation is enhanced during cell division, at which time vimentin filaments are significantly reorganized. Phosphorylation by PKN1 inhibits the formation of filaments. Phosphorylated at Ser-56 by CDK5 during neutrophil secretion in the cytoplasm. Phosphorylated by STK33. O-glycosylated during cytokinesis at sites identical or close to phosphorylation sites, this interferes with the phosphorylation status. S-nitrosylation is induced by interferon-gamma and oxidatively-modified low-densitity lipoprotein (LDL(ox)) possibly implicating the iNOS-S100A8/9 transnitrosylase complex.	
Cellular localization	Cytoplasm.	
Form	Vimentin is found in connective tissue and in the cytoskeleton.	
Images		




Flow Cytometry (Intracellular) - PE Anti-Vimentin antibody [EPR3776] - Cytoskeleton Marker (ab209446)

Overlay histogram showing HAP1 wildtype (green line) and HAP1-VIM knockout cells (red line) stained with ab209446. The cells were fixed with 80% methanol (5 min) and then permeabilized with 0.1% PBS-Triton X-100 for 15 min. The cells were then incubated in 1x PBS / 10% normal goat serum to block non-specific protein-protein interactions followed by the antibody (ab209446, 0.01 μ g/1x106 cells OR 1/5000 dilution) for 30 min at 22°C.

A rabbit IgG isotype control antibody (ab209478) was used at the same concentration and conditions as the primary antibody (HAP1 wildtype - black line, HAP1-VIM knockout - grey line). Unlabelled sample was also used as a control (this line is not shown for the purpose of simplicity).

Acquisition of >5,000 events were collected using a 50 mW Yellow/Green laser (561nm) and 586/15 bandpass filter.




Immunocytochemistry/ Immunofluorescence - PE Anti-Vimentin antibody [EPR3776] - Cytoskeleton Marker (ab209446)

ab209446 staining Vimentin in HeLa cells. The cells were fixed with 100% methanol (5min), permeabilized with 0.1% Triton X-100 for 5 minutes and then blocked with 1% BSA/10% normal goat serum/0.3M glycine in 0.1% PBS-Tween for 1h. The cells were then incubated overnight at +4°C with ab209446 at 1/500 dilution (Pseudocolored in green) and **ab195884**, Rat monoclonal to Tubulin (Alexa Fluor[®] 647), at 1/250 dilution (shown in red). Nuclear DNA was labelled with DAPI (shown in blue).

Image was taken with a confocal microscope (Leica-Microsystems, $TCS\ SP8$).

This product also gave a positive signal under the same testing conditions in HeLa cells fixed with 4% formaldehyde (10 min).




Flow Cytometry (Intracellular) - PE Anti-Vimentin antibody [EPR3776] - Cytoskeleton Marker (ab209446)

Overlay histogram showing HeLa cells stained with ab209446 (red line). The cells were fixed with 4% formaldehyde (10 min) and then permeabilized with 90% methanol (-20°C) for 30 min. The cells were then incubated in 1x PBS / 10% normal goat serum to block non-specific protein-protein interactions followed by the antibody (ab209446, 1/5000 dilution) for 30 min at 22°C.

Isotype control antibody (black line) was rabbit IgG (monoclonal) Phycorythrin (<u>ab209478</u>) used at the same concentration and conditions as the primary antibody. Unlabelled sample (blue line) was also used as a control.

Acquisition of >5,000 events were collected using a 20 mW Solid State Blue Laser (488nm) and 585/40 bandpass filter.



Cytoskeleton Marker (ab209446)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- · Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit https://www.abcam.com/abpromise or contact our technical team.

Terms and conditions

• Guarantee only valid for products bought direct from Abcam or one of our authorized distributors