# abcam

# Product datasheet

# Human FGF10 peptide ab91424

**Description** 

Product name Human FGF10 peptide

Purity > 70 % HPLC.

70 - 90% by HPLC

Animal free No

Nature Synthetic

**Species** Human

#### **Specifications**

Our Abpromise guarantee covers the use of ab91424 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Form Liquid

Additional notes - First try to dissolve a small amount of peptide in either water or buffer. The more charged

residues on a peptide, the more soluble it is in aqueous solutions.

- If the peptide doesn't dissolve try an organic solvent e.g. DMSO, then dilute using water or

buffer.

- Consider that any solvent used must be compatible with your assay. If a peptide does not

dissolve and you need to recover it, lyophilise to remove the solvent.

- Gentle warming and sonication can effectively aid peptide solubilisation. If the solution is

cloudy or has gelled the peptide may be in suspension rather than solubilised.

- Peptides containing cysteine are easily oxidised, so should be prepared in solution just prior

to use.

### Preparation and Storage

**Stability and Storage** Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw

cycles.

Information available upon request.

#### **General Info**

**Function** Could be a growth factor active in the process of wound healing. Acts as a mitogen in the lung.

1

May act in a manner similar to FGF-7.

#### Involvement in disease

Defects in FGF10 are the cause of autosomal dominant aplasia of lacrimal and salivary glands (ALSG) [MIM:180920]. ALSG has variable expressivity, and affected individuals may have aplasia or hypoplasia of the lacrimal, parotid, submandibular and sublingual glands and absence of the lacrimal puncta. The disorder is characterized by irritable eyes, recurrent eye infections, epiphora (constant tearing) and xerostomia (dryness of the mouth), which increases the risk of dental erosion, dental caries, periodontal disease and oral infections.

Defects in FGF10 are a cause of lacrimo-auriculo-dento-digital syndrome (LADDS)

[MIM:149730]; also known as Levy-Hollister syndrome. LADDS is a form of ectodermal dysplasia, a heterogeneous group of disorders due to abnormal development of two or more ectodermal structures. LADDS is an autosomal dominant syndrome characterized by aplastic/hypoplastic lacrimal and salivary glands and ducts, cup-shaped ears, hearing loss, hypodontia and enamel hypoplasia, and distal limb segments anomalies. In addition to these cardinal features, facial dysmorphism, malformations of the kidney and respiratory system and abnormal genitalia have been reported. Craniosynostosis and severe syndactyly are not observed.

Sequence similarities

Belongs to the heparin-binding growth factors family.

**Cellular localization** 

Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- · Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <a href="https://www.abcam.com/abpromise">https://www.abcam.com/abpromise</a> or contact our technical team.

## Terms and conditions

· Guarantee only valid for products bought direct from Abcam or one of our authorized distributors