

Product datasheet

Mineral Oil High sensitivity ab243855

Overview

Product name	Mineral Oil High sensitivity
General notes	<p>Ready to use as supplied. Pre-warm to 37°C prior to use.</p> <p>Although mineral oil is provided in a dropper bottle for convenience, we recommend using a repeater pipette for routine use. To apply oil using a repeater pipette, trim 3 – 4 mm off the tip at a 45° angle. Remove internal nozzle cap from the dropper bottle and slowly pick up the pre-warmed mineral oil.</p> <p>Store mineral oil at room temperature in the dark.</p> <p>ab243855 is a component of Mitochondrial Stress Test Complete Assay Kit (ab232857) and Extracellular Oxygen Consumption Assay (ab197243). It can be used in combination with Extracellular Oxygen Consumption Reagent (ab197242).</p>

Properties

Form	Liquid
Storage instructions	Shipped at Room Temperature. Store at Room Temperature. Store at +4°C. Store In the Dark.
Storage buffer	Constituents: 80% Butene Homopolymer, 5% Terphenyl

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors