abcam

Product datasheet

Silicon Imaging Chamber (1.75 mm) ab243309

* * * * * 1 Abreviews

Overview

Product name Silicon Imaging Chamber (1.75 mm)

Tested applications Suitable for: Histology

General notes After you have finished staining and clearing your tissues using Abcam's Tissue Clearing Kits and

reagents, you can use Silicon Imaging Chambers to mount your tissue for imaging. Silicon Imaging Chambers come in three different thicknesses and naturally adhere to glass microscope

slides. Simply place a cover slip on top and start imaging. Silicon Imaging Chamber are

chemically inert and can be used with any tissue clearing technique.

ab243309 - 1.75mm ab243311 - 3.5 mm ab243310 - 7mm

To learn more about tissue clearing, and tissue clearing methods, please review our guide to

tissue clearing.

Tissue clearing reagents and companion products are also available as individual products.

Properties

Form Solid

Storage instructions Shipped at Room Temperature. Store at Room Temperature.

Storage buffer Constituent: 100% Silicone oil

Applications

The Abpromise guarantee Our Abpromise guarantee covers the use of ab243309 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Histology		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- · Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- · We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit https://www.abcam.com/abpromise or contact our technical team.

Terms and conditions

• Guarantee only valid for products bought direct from Abcam or one of our authorized distributors