

**ab178078 –  
Mouse (FBV strain)  
tissue array of 22 organs,  
54 samples (1.5mm)**

**Instructions for Use**

Designed for IHC or ISH based protein or RNA tissue profiling in majority of Mouse organs.


This product is for research use only and is not intended for diagnostic use.

### Product Overview:

ab178078 is a normal tissue array of 22 major organs duplicated or triplicated from 3 FBV strain mice.

Case:	3
Core:	54
Format:	6 rows x 9 columns
Size:	1.5mm
Thickness:	4µm

All the tissues were fixed in 10% neutral buffered formalin for 24 hours and processed using identical SOPs. Sections were picked onto Superfrost Plus or Startfrost Adhesive slides. There may be >5% tissue core losses for a slide. But the core retention rate should be 90%.



### Notes for use:

Bake the slides at 60°C for 30 minutes before use. If antigen retrieval is needed, it is always a good idea to start with a protocol with weak to mild strength.

**For hints and tips on Immunohistochemistry as well as other techniques please visit [www.abcam.com/protocols](http://www.abcam.com/protocols).**

**Keep up to date on Abcam events and conferences on Facebook & Twitter.**

Array Position	Age	Sex	Anatomic Site	Nature	Pathology	Grade	Stage (TNM)
A01	3M	M	Bladder	Normal	Normal		
A02	3M	M	Bladder	Normal	Normal		
A03	3M	M	Brain	Normal	Normal		
A04	3M	F	Brain	Normal	Normal		
A05	3M	M	Brain	Normal	Normal		
A06	3M	M	Cerebellum	Normal	Normal		
A07	3M	F	Cerebellum	Normal	Normal		
A08	3M	M	Cerebellum	Normal	Normal		
A09	3M	F	Eye	Normal	Normal		
B01	3M	M	Esophagus	Normal	Normal		
B02	3M	F	Esophagus	Normal	Normal		
B03	3M	M	Stomach	Normal	Normal		
B04	3M	F	Stomach	Normal	Normal		
B05	3M	M	Stomach	Normal	Normal		
B06	3M	M	Intestine, small intestine	Normal	Normal		
B07	3M	F	Intestine, small intestine	Normal	Normal		
B08	3M	M	Intestine, colon	Normal	Normal		
B09	3M	F	Intestine, colon	Normal	Normal		
C01	3M	M	Heart	Normal	Normal		
C02	3M	F	Heart	Normal	Normal		
C03	3M	M	Heart	Normal	Normal		
C04	3M	M	Kidney	Normal	Normal		
C05	3M	F	Kidney	Normal	Normal		
C06	3M	M	Kidney	Normal	Normal		
C07	3M	M	Liver	Normal	Normal		
C08	3M	F	Liver	Normal	Normal		

Array Position	Age	Sex	Anatomic Site	Nature	Pathology	Grade	Stage (TNM)
C09	3M	M	Liver	Normal	Normal		
D01	3M	M	Lung	Normal	Normal		
D02	3M	F	Lung	Normal	Normal		
D03	3M	M	Lung	Normal	Normal		
D04	3M	F	Ovary	Normal	Normal		
D05	3M	F	Uterus	Normal	Normal		
D06	3M	F	Uterus	Normal	Normal		
D07	3M	M	Pancreas	Normal	Normal		
D08	3M	F	Pancreas	Normal	Normal		
D09	3M	M	Pancreas	Normal	Normal		
E01	3M	M	Prostate	Normal	Normal		
E02	3M	M	Prostate	Normal	Normal		
E03	3M	M	Salivary gland	Normal	Normal		
E04	3M	F	Salivary gland	Normal	Normal		
E05	3M	M	Salivary gland	Normal	Normal		
E06	3M	M	Skin	Normal	Normal		
E07	3M	F	Skin	Normal	Normal		
E08	3M	F	Skin	Normal	Normal		
E09	3M	M	Spleen	Normal	Normal		
F01	3M	F	Spleen	Normal	Normal		
F02	3M	M	Spleen	Normal	Normal		
F03	3M	M	Striated muscle	Normal	Normal		
F04	3M	M	Striated muscle	Normal	Normal		
F05	3M	F	Striated muscle	Normal	Normal		
F06	3M	M	Testis	Normal	Normal		
F07	3M	M	Testis	Normal	Normal		

Array Position	Age	Sex	Anatomic Site	Nature	Pathology	Grade	Stage (TNM)
F08	3M	M	Thymus	Normal	Normal		
F09	3M	M	Thymus	Normal	Normal		





**UK, EU and ROW**

Email: [technical@abcam.com](mailto:technical@abcam.com) | Tel: +44-(0)1223-696000

**Austria**

Email: [wissenschaftlicherdienst@abcam.com](mailto:wissenschaftlicherdienst@abcam.com) | Tel: 019-288-259

**France**

Email: [supportscientifique@abcam.com](mailto:supportscientifique@abcam.com) | Tel: 01-46-94-62-96

**Germany**

Email: [wissenschaftlicherdienst@abcam.com](mailto:wissenschaftlicherdienst@abcam.com) | Tel: 030-896-779-154

**Spain**

Email: [soportecientifico@abcam.com](mailto:soportecientifico@abcam.com) | Tel: 911-146-554

**Switzerland**

Email: [technical@abcam.com](mailto:technical@abcam.com)

Tel (Deutsch): 0435-016-424 | Tel (Français): 0615-000-530

**US and Latin America**

Email: [us.technical@abcam.com](mailto:us.technical@abcam.com) | Tel: 888-77-ABCAM (22226)

**Canada**

Email: [ca.technical@abcam.com](mailto:ca.technical@abcam.com) | Tel: 877-749-8807

**China and Asia Pacific**

Email: [hk.technical@abcam.com](mailto:hk.technical@abcam.com) | Tel: 108008523689 (中國聯通)

**Japan**

Email: [technical@abcam.co.jp](mailto:technical@abcam.co.jp) | Tel: +81-(0)3-6231-0940

[www.abcam.com](http://www.abcam.com) | [www.abcam.cn](http://www.abcam.cn) | [www.abcam.co.jp](http://www.abcam.co.jp)