

Product datasheet

Anti-ABCF3 antibody ab82478

Overview

**Product name** Anti-ABCF3 antibody  
**Description** Goat polyclonal to ABCF3

**i** This product is a **fast track antibody**. It has been affinity purified and shows high titre values against the immunizing peptide by ELISA.

[Read the terms of use »](#)

**Host species** Goat

**Species reactivity**

**Predicted to work with:** Mouse, Rat, Cow, Human, Xenopus laevis **!**

**Immunogen**

Synthetic peptide:

C-KESRLESSGKNKSYD

, corresponding to internal sequence amino acids 163-177 of human ABCF3 (NP\_060828.2).

[Run BLAST with](#)

[Run BLAST with](#)

Properties

**Form** Liquid  
**Storage instructions** Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.  
**Storage buffer** Preservative: 0.02% Sodium Azide  
 Constituents: 0.5% BSA, Tris saline, pH 7.3  
**Purity** Immunogen affinity purified  
**Purification notes** Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.  
**Clonality** Polyclonal  
**Isotype** IgG

Applications

**Application notes** Peptide ELISA: antibody detection limit dilution 1:128,000.  
 WB: Preliminary experiments in human liver lysates and in lysates of cell lines DAUDI and K562 gave no specific signal but low background (at antibody concentration up to 1 µg/ml). We would

welcome your feedback if you have tested this antibody in more suitable positive controls.

## Target

---

### Relevance

ABCF3 belongs to the ABC transporter superfamily. ABCF family. EF3 subfamily. Contains 2 ABC transporter domains. There are 2 named isoforms produced by alternative splicing.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors