

Product datasheet

Anti-ABCF3 antibody ab82478

Overview

Product name Anti-ABCF3 antibody
Description Goat polyclonal to ABCF3

i This product is a **fast track antibody**. It has been affinity purified and shows high titre values against the immunizing peptide by ELISA.

[Read the terms of use »](#)

Species reactivity

Predicted to work with: Mouse, Rat, Cow, Human, Xenopus laevis 

Immunogen

Synthetic peptide:

C-KESRLESSGKNKSYD

, corresponding to internal sequence amino acids 163-177 of human ABCF3 (NP_060828.2).

 [Run BLAST with](#)

 [Run BLAST with](#)

Properties

Form Liquid
Storage instructions Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
Storage buffer Preservative: 0.02% Sodium Azide
 Constituents: 0.5% BSA, Tris saline, pH 7.3
Purity Immunogen affinity purified
Purification notes Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Clonality Polyclonal
Isotype IgG

Applications

Application notes Peptide ELISA: antibody detection limit dilution 1:128,000.

WB: Preliminary experiments in human liver lysates and in lysates of cell lines DAUDI and K562 gave no specific signal but low background (at antibody concentration up to 1 µg/ml). We would welcome your feedback if you have tested this antibody in more suitable positive controls.

Target

Relevance

ABCF3 belongs to the ABC transporter superfamily. ABCF family. EF3 subfamily. Contains 2 ABC transporter domains. There are 2 named isoforms produced by alternative splicing.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors