

Product datasheet

Anti-APG8A antibody ab98830

1 References

Overview

Product name	Anti-APG8A antibody
Description	Rabbit polyclonal to APG8A
Host species	Rabbit
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Arabidopsis thaliana
Immunogen	Full length protein produced in E.coli.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Preservative: None Constituents: Whole serum
Purity	Whole antiserum
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab98830** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/1000. Predicted molecular weight: 14 kDa.

Target

Relevance	APG8A is involved in cytoplasm to vacuole transport (Cvt) vesicles and autophagosomes formation. It may mediate the delivery of the vesicles and autophagosomes to the vacuole via the
------------------	--

microtubule cytoskeleton.

Cellular localization

Cytoplasmic vesicle › cvt vesicle membrane; Lipid-anchor. Cytoplasmic vesicle › autophagosome membrane; Lipid-anchor. Vacuole membrane; Lipid-anchor. Note= Membrane-associated through a lipid anchor.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors