

Product datasheet

Anti-APG8A antibody ab98830

1 References

Overview

Product name	Anti-APG8A antibody
Description	Rabbit polyclonal to APG8A
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Arabidopsis thaliana
Immunogen	Full length protein produced in E.coli.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Preservative: None Constituents: Whole serum
Purity	Whole antiserum
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab98830** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/1000. Predicted molecular weight: 14 kDa.

Target

Relevance	APG8A is involved in cytoplasm to vacuole transport (Cvt) vesicles and autophagosomes formation. It may mediate the delivery of the vesicles and autophagosomes to the vacuole via the microtubule cytoskeleton.
------------------	--

Cellular localization

Cytoplasmic vesicle › cvt vesicle membrane; Lipid-anchor. Cytoplasmic vesicle › autophagosome membrane; Lipid-anchor. Vacuole membrane; Lipid-anchor. Note= Membrane-associated through a lipid anchor.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors