

Product datasheet

Anti-C11orf17 antibody ab69697

1 Image

Overview

Product name	Anti-C11orf17 antibody
Description	Mouse polyclonal to C11orf17
Host species	Mouse
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human
Immunogen	Full length protein (Human) C11orf17 (NP_065693).
Positive control	C11orf17 transfected 293T cell lysate

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Preservative: None Constituents: 1X PBS, pH 7.2
Purity	Protein G purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab69697** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


Application	Abreviews	Notes
WB		1/500 - 1/1000. Detects a band of approximately 26 kDa (predicted molecular weight: 23 kDa).

Target

Relevance C11orf17 is a nuclear protein and known to bind to the amino terminus (residues 1-39) of the C subunit of PKA. It is a binding partner of NF-kappaB p65 subunit, and enhances NF-kappaB-mediated gene expression.

Cellular localization Nuclear

Images



Western blot - C11orf17 antibody (ab69697)

All lanes : Anti-C11orf17 antibody (ab69697)
at 1/500 dilution

Lane 1 : C11orf17 transfected 293T cell
lysate

Lane 2 : Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-Mouse IgG (H&L)-HRP
Conjugate at 1/2500 dilution

Predicted band size: 23 kDa

Observed band size: 26 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors