

## Product datasheet

# Anti-CCNO antibody ab47681

★★★★☆ 1 Abreviews 1 Image

### Overview

<b>Product name</b>	Anti-CCNO antibody
<b>Description</b>	Rabbit polyclonal to CCNO
<b>Specificity</b>	At least two isoforms of Cyclin O are known to exist; this antibody will recognize both isoforms.
<b>Tested applications</b>	<b>Suitable for:</b> WB
<b>Species reactivity</b>	<b>Reacts with:</b> Mouse, Rat, Human
<b>Immunogen</b>	A 13 amino acid peptide from near the center of Human Cyclin O (Genbank accession No. NP_066970).
<b>Positive control</b>	Human bladder tissue lysate.

### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C.
<b>Storage buffer</b>	Preservative: 0.02% Sodium Azide Constituents: PBS
<b>Purity</b>	Immunogen affinity purified
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

### Applications

Our [Abpromise guarantee](#) covers the use of **ab47681** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use a concentration of 1 - 2 µg/ml. Detects a band of approximately 40 kDa (predicted molecular weight: 38 kDa).

### Target


**Sequence similarities**

Belongs to the cyclin family.

**Developmental stage**

Maximum levels during G(1) phase. Levels decrease through S and G(2) phases.

---

**Images**

**Predicted band size : 38 kDa**

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors