

## Product datasheet

# Anti-CD30 antibody [mCD30.1] (Phycoerythrin) ab95650

### Overview

<b>Product name</b>	Anti-CD30 antibody [mCD30.1] (Phycoerythrin)
<b>Description</b>	Armenian Hamster monoclonal [mCD30.1] to CD30 (Phycoerythrin)
<b>Host species</b>	Armenian hamster
<b>Conjugation</b>	Phycoerythrin. Ex: 488nm, Em: 575nm
<b>Tested applications</b>	<b>Suitable for:</b> Flow Cyt
<b>Species reactivity</b>	<b>Reacts with:</b> Mouse
<b>Immunogen</b>	Fusion protein of Mouse CD30-IgG1
<b>Positive control</b>	Activated Mouse splenocyte suspensions

### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C.
<b>Storage buffer</b>	Preservative: 0.09% Sodium Azide Constituents: PBS, 150mM Sodium chloride, pH 7.2
<b>Purity</b>	Protein G purified
<b>Clonality</b>	Monoclonal
<b>Clone number</b>	mCD30.1
<b>Isotype</b>	IgG

### Applications

Our [Abpromise guarantee](#) covers the use of **ab95650** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt		Use 0.5µg for 10 <sup>5-8</sup> cells.

### Target

<b>Function</b>	Receptor for TNFSF8/CD30L. May play a role in the regulation of cellular growth and transformation of activated lymphoblasts. Regulates gene expression through activation of NF-kappa-B.
<b>Sequence similarities</b>	Contains 6 TNFR-Cys repeats.
<b>Post-translational modifications</b>	Phosphorylated on serine and tyrosine residues.
<b>Cellular localization</b>	Cytoplasm and Cell membrane.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors