

Product datasheet

Anti-CD5 antibody [M28623] (FITC) ab52310

1 References

Overview

Product name	Anti-CD5 antibody [M28623] (FITC)
Description	Mouse monoclonal [M28623] to CD5 (FITC)
Host species	Mouse
Conjugation	FITC. Ex: 493nm, Em: 528nm
Specificity	This antibody reacts with CD5.
Tested applications	Suitable for: Flow Cyt, ICC/IF
Species reactivity	Reacts with: Human
Immunogen	The details of the immunogen for this antibody are not available.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: 0.08% Sodium Azide Constituents: PBS, 0.2% protein carrier
Clonality	Monoclonal
Clone number	M28623
Isotype	IgG2a
Light chain type	kappa

Applications

Our [Abpromise guarantee](#) covers the use of **ab52310** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt		
ICC/IF		

Application notes

Flow Cyt: Use 10µl for 10⁶ cells.

Flow Cyt: Use 10µl for 100µl whole blood.

ICC/IF: Use at an assay dependent dilution.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Target

Function

May act as a receptor in regulating T-cell proliferation. CD5 interacts with CD72/LYB-2.

Sequence similarities

Contains 3 SRCR domains.

Cellular localization

Membrane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors