

Product datasheet

COLO205 whole cell lysate ab95356

Overview

| | |
|----------------------------|---|
| Product name | COLO205 whole cell lysate |
| General notes | <p>Growth Properties: mixed, adherent and suspension</p> <p>Morphology: epithelial</p> <p>Organ: colon</p> <p>Tumor Stage: Dukes' type D</p> <p>Disease: colorectal adenocarcinoma</p> <p>Derived from metastatic site: ascites</p> <p>Cellular Products: carcinoembryonic antigen (CEA) 1.5 to 4.1 ng/10 exp6 cells/10 days; keratin; interleukin 10 (IL-10, interleukin-10)</p> <p>Tumorigenic: Yes</p> <p>This is a recommended negative control for PD-L1 RabMAb (clone 28-8), Ab205921</p> |
| Tested applications | Suitable for: WB, SDS-PAGE |

Properties

| | |
|-----------------------------|---|
| Form | Lyophilised: Add 100 µl distilled water to the sample. Final concentration is 2 mg/ml. |
| Storage instructions | Shipped at 4°C. After reconstitution store at -20°C. Avoid freeze / thaw cycles. |
| Storage buffer | pH: 6.80 Constituents: 0.1% Bromophenol blue, 2% SDS, 10% Glycerol, 1.54% DTT, 0.79% Tris HCl |
| Purity | Whole Cell Lysate |
| Lysate notes | This line was isolated in 1975 by T.U. Semple, et al. from ascitic fluid of a 70-year-old Caucasian male with carcinoma of the colon. The patient had been treated with 5-fluorouracil for 4-6 weeks before removal of the fluid specimen. The cells are CSAp negative (CSAp-). The line was derived from tissue from the same patient as COLO 201. The cells are positive for keratin by immunoperoxidase staining. COLO 205 cells express a 36000 dalton cell surface glycoprotein related to the GA733-2 tumor associated antigen. |

Applications

Our [Abpromise guarantee](#) covers the use of **ab95356** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| Application | Abreviews | Notes |
|-------------|-----------|--|
| WB | | Use at an assay dependent concentration. |
| SDS-PAGE | | Use at an assay dependent concentration. |

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors