

## Product datasheet

# Anti-Creatinine antibody ab30719

### 1 References

#### Overview

<b>Product name</b>	Anti-Creatinine antibody
<b>Description</b>	Sheep polyclonal to Creatinine
<b>Host species</b>	Sheep
<b>Specificity</b>	Recognizes free creatinine.
<b>Tested applications</b>	<b>Suitable for:</b> Competitive ELISA
<b>Immunogen</b>	Creatinine(N)-BSA

#### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
<b>Storage buffer</b>	Preservative: 0.09% Sodium Azide Constituents: 20mM PBS, 150mM Sodium chloride, pH 7.2
<b>Purity</b>	IgG fraction
<b>Purification notes</b>	Prepared by caprylic acid and ammonium sulphate precipitation techniques.
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

#### Applications

Our [Abpromise guarantee](#) covers the use of **ab30719** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Competitive ELISA		1/6000.

#### Target

**Relevance** Creatinine is a waste product from protein in the diet and from the muscles of the body.

Creatinine is removed from the body by the kidneys; as kidney disease progresses, the level of creatinine in the blood increases.

**Cellular localization**

Secreted

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

**Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors