

Product datasheet

D Serine small molecule ab102605

Overview

Product name D Serine small molecule

Description

Nature Recombinant

Amino Acid Sequence

Specifications

Our [Abpromise guarantee](#) covers the use of **ab102605** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Form Lyophilised

Preparation and Storage

Stability and Storage Shipped at 4°C.

Reconstitution Lyophilized antigen are stable at least 2 years. For reconstitution, use pure water with 0.1% merthiolate. This solution is stable at +4°C for 2 months.

General Info

Cellular localization Ubiquitous

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors