

## Product datasheet

# D Serine small molecule ab102605

### Overview

---

**Product name** D Serine small molecule

### Description

---

**Nature** Recombinant

### Amino Acid Sequence

### Specifications

---

Our [Abpromise guarantee](#) covers the use of **ab102605** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**Form** Lyophilised

### Preparation and Storage

---

**Stability and Storage** Shipped at 4°C.

**Reconstitution** Lyophilized antigen are stable at least 2 years. For reconstitution, use pure water with 0.1% merthiolate. This solution is stable at +4°C for 2 months.

### General Info

---

**Cellular localization** Ubiquitous

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors