

Product datasheet

Anti-DR3 antibody ab8411

1 Image

Overview

Product name	Anti-DR3 antibody
Description	Rabbit polyclonal to DR3
Specificity	Does not recognise other death receptors.
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Mouse, Rat, Human
Immunogen	Synthetic peptide, corresponding to C terminal amino acids 398-417 of Human DR3 (Peptide available as ab8443 .)
Positive control	Jurkat total cell lysate A 59 kDa band should be detected.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	PBS with 0.02% sodium azide
Purity	IgG fraction
Clonality	Polyclonal
Isotype	IgG
Light chain type	unknown

Applications

Our [Abpromise guarantee](#) covers the use of **ab8411** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/500 - 1/1000. Detects a band of approximately 59 kDa. Can be blocked with Human DR3 peptide (ab8443) .

Target

Function	Receptor for TNFSF12/APO3L/TWEAK. Interacts directly with the adapter TRADD. Mediates
-----------------	---

activation of NF-kappa-B and induces apoptosis. May play a role in regulating lymphocyte homeostasis.

Tissue specificity

Abundantly expressed in thymocytes and lymphocytes. Detected in lymphocyte-rich tissues such as thymus, colon, intestine, and spleen. Also found in the prostate.

Sequence similarities

Contains 1 death domain.
Contains 4 TNFR-Cys repeats.


Post-translational modifications

Glycosylated.

Cellular localization

Secreted and Cell membrane.

Anti-DR3 antibody images



Anti-DR3 antibody (ab8411) at 1/500 dilution
+ Jurkat total cell lysate

Observed band size : 59 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors