

Product datasheet

Donkey Anti-Mouse IgG H&L (Alkaline Phosphatase) preadsorbed ab7062

Overview

Product name	Donkey Anti-Mouse IgG H&L (Alkaline Phosphatase) preadsorbed	
Target species	Mouse	
Tested applications	Suitable for: WB, ICC/IF, Dot blot, ELISA, IHC-P, IHC-Fr, Immunomicroscopy	
Minimal cross-reactivity	Chicken, Cow, Goat, Guinea pig, Hamster, Horse, Human, Rabbit, Rat, Sheep	more details
Immunogen	Full length native mouse IgG (purified).	
Conjugation	Alkaline Phosphatase	

Properties

Form	Liquid	
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.	
Storage buffer	Preservative: 0.1% Sodium Azide Constituents: 50% Glycerol, 10mg/ml BSA, 0.15M Sodium chloride, 0.05M Tris HCl, 1mM Magnesium chloride, 0.1mM Zinc chloride, pH 8.0	
Purity	IgG fraction	
Purification notes	This product was prepared from monospecific antiserum by immunoaffinity chromatography using mouse IgG coupled to agarose beads.	
Conjugation notes	Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons)	
Clonality	Polyclonal	

Applications

Our [Abpromise guarantee](#) covers the use of **ab7062** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		
ICC/IF		

Application	Abreviews	Notes
Dot blot		
ELISA		
IHC-P		
IHC-Fr		
Immunomicroscopy		

Application notes

Dot: Use at an assay dependant dilution.

ELISA: 1/500 - 1/2000.

IHC-P: Use at an assay dependant dilution.

IHC-Fr: Use at an assay dependant dilution.

IM: Use at an assay dependant dilution.

WB: Use at an assay dependant dilution.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors