

Product datasheet

Donkey Anti-Rabbit IgG H&L (DyLight® 594) preadsorbed ab96921

[4 References](#) [2 Images](#)

Overview

Product name	Donkey Anti-Rabbit IgG H&L (DyLight® 594) preadsorbed
Description	Donkey polyclonal Secondary Antibody to Rabbit IgG - H&L (DyLight® 594), pre-adsorbed
Target species	Rabbit
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Rabbit IgG and with light chains common to other Rabbit immunoglobulins. No antibody was detected against non immunoglobulin serum proteins. Less than 0.1% cross reactivity to bovine, chicken, goat, horse, human, mouse, pig and rat was detected.
Tested applications	Suitable for: Flow Cyt, ICC/IF, IHC-P
Minimal cross-reactivity	Chicken, Cow, Goat, Horse, Human, Mouse, Pig, Rat more details
Conjugation	DyLight® 594. Ex: 593nm, Em: 618nm

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: 0.09% Sodium Azide Constituents: 0.2% BSA, PBS
Purity	Immunogen affinity purified
Purification notes	Antiserum was cross adsorbed using bovine, chicken, goat, human, mouse, pig and rat immunosorbents to remove cross reactive antibodies. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to DyLight® 594.
Clonality	Polyclonal
Isotype	IgG

Applications


Our [Abpromise guarantee](#) covers the use of **ab96921** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
-------------	-----------	-------

Flow Cyt	1/50 - 1/200.
ICC/IF	1/50 - 1/500.
IHC-P	1/50 - 1/500.


Donkey Anti-Rabbit IgG H&L (DyLight® 594) preadsorbed images



E14 mouse brain sample was stained with primary antibody Rabbit anti-Clim1. The secondary antibody used was ab96921 (1:200 dilution, incubate for 1 hour at room temperature). The blue background was from DAPI.

This image is courtesy of an anonymous Abreview

Immunohistochemistry (Frozen sections) - Donkey polyclonal Secondary Antibody to Rabbit IgG - H&L (DyLight® 594), pre-adsorbed (ab96921)



Donkey Anti-Rabbit IgG H&L (DyLight® 594) preadsorbed (ab96921)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors