

## Product datasheet

# Donkey Anti-Rabbit IgG H&L (DyLight® 594) preadsorbed ab96921

[6 References](#)   [2 Images](#)

### Overview

---

<b>Product name</b>	Donkey Anti-Rabbit IgG H&L (DyLight® 594) preadsorbed
<b>Description</b>	Donkey polyclonal Secondary Antibody to Rabbit IgG - H&L (DyLight® 594), pre-adsorbed
<b>Target species</b>	Rabbit
<b>Specificity</b>	By immunoelectrophoresis and ELISA this antibody reacts specifically with Rabbit IgG and with light chains common to other Rabbit immunoglobulins. No antibody was detected against non immunoglobulin serum proteins. Less than 0.1% cross reactivity to bovine, chicken, goat, horse, human, mouse, pig and rat was detected.
<b>Tested applications</b>	<b>Suitable for:</b> Flow Cyt, ICC/IF, IHC-P
<b>Minimal cross-reactivity</b>	Chicken, Cow, Goat, Horse, Human, Mouse, Pig, Rat <a href="#">more details</a>
<b>Conjugation</b>	DyLight® 594. Ex: 593nm, Em: 618nm

### Properties

---

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C.
<b>Storage buffer</b>	Preservative: 0.09% Sodium Azide Constituents: 0.2% BSA, PBS
<b>Purity</b>	Immunogen affinity purified
<b>Purification notes</b>	Antiserum was cross adsorbed using bovine, chicken, goat, human, mouse, pig and rat immunosorbents to remove cross reactive antibodies. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to DyLight® 594.
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

### Applications


---

Our [Abpromise guarantee](#) covers the use of **ab96921** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt		1/50 - 1/200.
ICC/IF		1/50 - 1/500.
IHC-P		1/50 - 1/500.


## Images



E14 mouse brain sample was stained with primary antibody Rabbit anti-Clim1. The secondary antibody used was ab96921 (1:200 dilution, incubate for 1 hour at room temperature). The blue background was from DAPI.

*This image is courtesy of an anonymous Abreview*

Immunohistochemistry (Frozen sections) - Donkey polyclonal Secondary Antibody to Rabbit IgG - H&L (DyLight® 594), pre-adsorbed (ab96921)



Donkey Anti-Rabbit IgG H&L (DyLight® 594)  
preadsorbed (ab96921)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors