

Product datasheet

Anti-FTHL17 antibody ab105955

1 Image

Overview

Product name	Anti-FTHL17 antibody
Description	Rabbit polyclonal to FTHL17
Host species	Rabbit
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human
Immunogen	Synthetic peptide, corresponding to a region within N-terminal amino acids 38-87 (MAFYFNRDDV ALENFFRYFL RLSDDKMEHA QKLMRLQNL RGGHICLHDIR) of Human FTHL17 (NP_114100). Run BLAST with ExPASy Run BLAST with NCBI
Positive control	721_B Cell Lysate

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Storage buffer	Preservative: None Constituents: 2% Sucrose, PBS
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab105955** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 21 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

Target


Tissue specificity

Testis specific. Also expressed in several cancers.

Sequence similarities

Belongs to the ferritin family.
Contains 1 ferritin-like diiron domain.

Images



Anti-FTHL17 antibody (ab105955) at 1 µg/ml
+ 721_B Cell Lysate at 10 µg

Predicted band size: 21 kDa

12% SDS PAGE

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors