


Product datasheet

Anti-Gm13178 antibody ab113822

1 Image

Overview

Product name	Anti-Gm13178 antibody
Description	Rabbit polyclonal to Gm13178
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Mouse Predicted to work with: Rat, Rabbit, Horse, Chicken, Guinea pig, Cow, Dog, Pig 
Immunogen	Synthetic peptide from a region within: TFLVSCEHDV LRDDALLYKK RLEDQGVVPS WYHAEDGFHG CISLFDKQPF, corresponding to C terminal amino acids 340-389 of Mouse Gm13178 (NP_001079005). Run BLAST with ExPASy Run BLAST with NCBI
Positive control	Mouse kidney lysate

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Constituents: 97% PBS, 2% Sucrose
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG


Applications

Our [Abpromise guarantee](#) covers the use of **ab113822** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 46 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

Images



Anti-Gm13178 antibody (ab113822) at 1
 $\mu\text{g/ml}$ + Mouse kidney lysate at 10 μg

Predicted band size : 46 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors