

Product datasheet

Goat Anti-Bird IgY H&L (FITC) ab112772

Overview

Product name	Goat Anti-Bird IgY H&L (FITC)
Description	Goat polyclonal Secondary Antibody to Bird IgY - H&L (FITC)
Target species	Bird
Specificity	By immunoelectrophoresis and ELISA ab112772 reacts specifically with Duck, Sparrow, Dove and Chicken IgG. Turkey, Ostrich, Emu, Quail, Macaw and Cockateil IgG are also detected. No antibody was detected against non-immunoglobulin serum proteins. ab112772 may cross react with IgG from other species.
Tested applications	Suitable for: IHC-P, ICC/IF
Conjugation	FITC. Ex: 493nm, Em: 528nm

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: 0.09% Sodium azide Constituents: 0.2% BSA, 99% PBS
Purity	Immunogen affinity purified
Purification notes	Goats were immunized with purified IgG from Dove, Duck, Sparrow and Chicken serum. ab112772 was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to fluorescein isothiocyanate (FITC).
Conjugation notes	F/P ratio is 4.4.
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab112772** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		1/200 - 1/2000.

Application	Abreviews	Notes
ICC/IF		1/200 - 1/2000.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors