

Product datasheet

Goat F(ab')₂ Anti-Mouse IgG - (Fab')₂ (DyLight® 550), pre-adsorbed ab98758

1 Image

Overview

Product name	Goat F(ab') ₂ Anti-Mouse IgG - (Fab') ₂ (DyLight® 550), pre-adsorbed
Host species	Goat
Target species	Mouse
Tested applications	Suitable for: Flow Cyt, ICC/IF, IHC-P
Minimal cross-reactivity	Human, Rat more details
Conjugation	DyLight® 550. Ex: 562nm, Em: 576nm

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: 0.09% Sodium Azide Constituents: 0.2% BSA, PBS
Purity	Immunogen affinity purified
Purification notes	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads. F(ab') ₂ fragment were generated using a pepsin digestion. Fc fragments and whole IgG molecules have been removed. Fragments were conjugated to DyLight® 550.
Clonality	Polyclonal
Isotype	IgG
General notes	DyLight® 550 replaces DyLight® 549.

Applications


Our [Abpromise guarantee](#) covers the use of **ab98758** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt		1/50 - 1/200.

Application	Abreviews	Notes
ICC/IF		1/50 - 1/500.
IHC-P		1/50 - 1/500.

Images



Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.

DyLight®-Goat F(ab')₂ polyclonal Secondary Antibody to Mouse IgG - (Fab')₂ (DyLight® 550), pre-adsorbed(ab98758)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors