

## Product datasheet

# Goat F(ab')<sub>2</sub> Anti-Mouse IgG H&L (SureLight® PE) ab99914

### Overview

|                            |   |
|----------------------------|---|
| <b>Product name</b>        | Goat F(ab') <sub>2</sub> Anti-Mouse IgG H&L (SureLight® PE) |
| <b>Host species</b>        | Goat  |
| <b>Target species</b>      | Mouse   |
| <b>Tested applications</b> | <b>Suitable for:</b> Flow Cyt                               |
| <b>Conjugation</b>         | SureLight® PE   |

### Properties

|                             |   |
|-----------------------------|---|
| <b>Form</b>                 | Lyophilised: Reconstitute to 0.5 ml with distilled deionized water, vortex gently and allow vial to sit on ice for 20 minutes.  |
| <b>Storage instructions</b> | Shipped at 4°C. Store at +4°C.  |
| <b>Storage buffer</b>       | Preservative: 0.05% Sodium Azide<br>Constituents: 0.1% BSA, 50mM Sucrose, 150mM Sodium chloride, 100mM Sodium phosphate, pH 7.4 |
| <b>Purity</b>               | IgG fraction  |
| <b>Clonality</b>            | Polyclonal  |
| <b>Isotype</b>              | IgG   |
| <b>General notes</b>        | Lyophilized material is stable for up to one year. After reconstitution should be used within 3 months.                         |

### Applications

Our [Abpromise guarantee](#) covers the use of **ab99914** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| Application | Abreviews | Notes                           |
|-------------|-----------|---------------------------------|
| Flow Cyt    |           | Use a concentration of 3 µg/ml. |

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors