

Product datasheet

Goat F(ab')₂ Anti-Mouse IgG H&L (SureLight® PE) ab99914

Overview

Product name	Goat F(ab') ₂ Anti-Mouse IgG H&L (SureLight® PE)
Host species	Goat
Target species	Mouse
Tested applications	Suitable for: Flow Cyt
Conjugation	SureLight® PE

Properties

Form	Lyophilised: Reconstitute to 0.5 ml with distilled deionized water, vortex gently and allow vial to sit on ice for 20 minutes.
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: 0.05% Sodium Azide Constituents: 0.1% BSA, 50mM Sucrose, 150mM Sodium chloride, 100mM Sodium phosphate, pH 7.4
Purity	IgG fraction
Clonality	Polyclonal
Isotype	IgG
General notes	Lyophilized material is stable for up to one year. After reconstitution should be used within 3 months.

Applications

Our [Abpromise guarantee](#) covers the use of **ab99914** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt		Use a concentration of 3 µg/ml.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors