

## Product datasheet

# Goat Anti-Mouse IgG1 heavy chain ab97236

### Overview

<b>Product name</b>	Goat Anti-Mouse IgG1 heavy chain
<b>Host species</b>	Goat
<b>Target species</b>	Mouse
<b>Specificity</b>	By immunoelectrophoresis and ELISA this antibody reacts specifically with mouse IgG1. Cross reactivity with other immunoglobulins is less than 2%.
<b>Tested applications</b>	<b>Suitable for:</b> ELISA, WB, IHC-P, ICC

### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	Preservative: 0.09% Sodium Azide Constituents: PBS
<b>Purity</b>	Immunogen affinity purified
<b>Purification notes</b>	Antiserum was solid phase adsorbed to ensure subclass specificity. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads.
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

### Applications

Our [Abpromise guarantee](#) covers the use of **ab97236** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		1/100 - 1/30000. (Coating ELISA: 1/100 - 1/500, Primary ELISA: 1/1000-1/30000)
WB		1/1000 - 1/30000. (Chemiluminescent: 1/1000 - 1/30000, Colorimetric: 1/1000 - 1/10000)
IHC-P		1/200 - 1/2000.

Application	Abreviews	Notes
ICC		Use at an assay dependent dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors