

Product datasheet

Goat Anti-Mouse IgG3 heavy chain (HRP) preadsorbed ab98708

Overview

Product name	Goat Anti-Mouse IgG3 heavy chain (HRP) preadsorbed	
Description	Goat polyclonal Secondary Antibody to Mouse IgG3 - heavy chain (HRP), pre-adsorbed	
Host species	Goat	
Target species	Mouse	
Tested applications	Suitable for: WB, IHC-P, ELISA, ICC	
Minimal cross-reactivity	Human, Rat	more details
Conjugation	HRP	

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: 0.1% Proclin Constituents: 0.2% BSA, PBS
Purity	Immunogen affinity purified
Purification notes	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Horse Radish Peroxidase.
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab98708** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/2000 - 1/25000. (Colorimetric: 1/2000 - 1/20000; Chemiluminescent: 1/5000 - 1/25000).
IHC-P		1/200 - 1/2000.

Application	Abreviews	Notes
ELISA		1/10000 - 1/100000. (as Primary).
ICC		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors