

Product datasheet

Goat Anti-Mouse IgM mu chain (6nm Gold) preadsorbed ab105292

Overview

Product name	Goat Anti-Mouse IgM mu chain (6nm Gold) preadsorbed
Target species	Mouse
Specificity	Recognizes rat IgM. Exhibits minimal crossreactivity to bovine, human and horse serum proteins.
Tested applications	Suitable for: Electron Microscopy, WB, ELISA, Other
Minimal cross-reactivity	Cow, Horse, Human more details
Conjugation	Gold 6nm

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: 0.05% Sodium Azide Constituents: 15mg/ml BSA, 0.15M Sodium chloride, 10mM Sodium phosphate, 10mM Sodium borate, pH 8.5
Purity	Immunogen affinity purified
Conjugation notes	6nm Gold
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab105292** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Electron Microscopy		1/10 - 1/20.
WB		Use at an assay dependent dilution.

Application	Abreviews	Notes
ELISA		Use at an assay dependent dilution.
Other		Use at an assay dependent dilution. Light microscopy

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors