

Product datasheet

Anti-Hsp70 antibody ab90554

1 Image

Overview

Product name	Anti-Hsp70 antibody
Description	Rabbit polyclonal to Hsp70
Host species	Rabbit
Tested applications	Suitable for: WB, IP
Species reactivity	Reacts with: Mouse, Rat, Sheep, Guinea pig, Hamster, Cow, Dog, Human, Pig, Carp, Fish, Monkey
Immunogen	Recombinant Human Hsp70 protein
Positive control	HeLa Cell Lysate and HeLa (Heat Shocked) Cell Lysate

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Preservative: 0.09% Sodium Azide Constituents: 50% Glycerol, PBS, pH 7.2
Purity	Protein A purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab90554** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/1000. Predicted molecular weight: 70 kDa.
IP		1/100.

Target

Function

In cooperation with other chaperones, Hsp70s stabilize preexistent proteins against aggregation and mediate the folding of newly translated polypeptides in the cytosol as well as within organelles. These chaperones participate in all these processes through their ability to recognize nonnative conformations of other proteins. They bind extended peptide segments with a net hydrophobic character exposed by polypeptides during translation and membrane translocation, or following stress-induced damage. In case of rotavirus A infection, serves as a post-attachment receptor for the virus to facilitate entry into the cell.

Tissue specificity

HSPA1B is testis-specific.


Sequence similarities

Belongs to the heat shock protein 70 family.

Cellular localization

Cytoplasm. Localized in cytoplasmic mRNP granules containing untranslated mRNAs.

Images



All lanes : Anti-Hsp70 antibody (ab90554) at 1/1000 dilution

Lane 1 : Molecular weight marker

Lane 2 : HeLa cell lysate

Lane 3 : HeLa cell lysate (Heat Shocked)

Predicted band size: 70 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors