

Product datasheet

Human ERp72 peptide ab39767

Overview

Product name Human ERp72 peptide

Description

Nature Synthetic
Amino Acid Sequence
Species Human
Sequence KFIEEHATKLSRTKEEL
Amino acids 629 to 645

Specifications

Our [Abpromise guarantee](#) covers the use of **ab39767** in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications Neutralising
 Blocking
Purity > 90 % SDS-PAGE.
Form Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

General Info

Sequence similarities Belongs to the protein disulfide isomerase family.
 Contains 3 thioredoxin domains.
Cellular localization Endoplasmic reticulum lumen. Melanosome. Identified by mass spectrometry in melanosome fractions from stage I to stage IV.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors