

## Product datasheet

# Human PRPF31 peptide **ab22909**

### Overview

---

**Product name** Human PRPF31 peptide

### Description

---

**Nature** Synthetic

#### Amino Acid Sequence

**Species** Human

**Sequence** C-LKVKGEKSGLMST

**Amino acids** 487 to 499

### Specifications

---

Our [Abpromise guarantee](#) covers the use of **ab22909** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**Applications** Blocking - Blocking peptide for Anti-PRPF31 antibody ([ab2467](#))

**Form** Liquid

### Preparation and Storage

---

**Stability and Storage** Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

### General Info

---

**Function** Involved in pre-mRNA splicing. Required for U4/U6.U5 tri-snRNP formation.

**Tissue specificity** Ubiquitously expressed.

**Involvement in disease** Defects in PRPF31 are the cause of retinitis pigmentosa type 11 (RP11) [MIM:600138]. RP leads to degeneration of retinal photoreceptor cells. Patients typically have night vision blindness and loss of midperipheral visual field. As their condition progresses, they lose their far peripheral visual field and eventually central vision as well. RP11 inheritance is autosomal dominant.

**Sequence similarities** Belongs to the PRP31 family.  
Contains 1 Nop domain.

<b>Domain</b>	Interacts with the snRNP via the Nop domain. The coiled coil domain is formed by two non-contiguous helices.
<b>Cellular localization</b>	Nucleus speckle. Nucleus > Cajal body. Predominantly found in speckles and in Cajal bodies.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

---

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

---

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors