

Product datasheet

Human RNF39 peptide ab45498

Overview

Product name Human RNF39 peptide

Description

Nature Synthetic

Amino Acid Sequence

Species Human

Sequence CDPRAPLRIVPAES

Amino acids 407 to 420

Specifications

Our [Abpromise guarantee](#) covers the use of **ab45498** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications Blocking - Blocking peptide for Anti-RNF39 antibody ([ab37898](#))

Form Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.

General Info

Function May play a role in prolonged long term-potential (LTP) maintenance.

Tissue specificity Expressed in testis.

Sequence similarities Contains 1 B30.2/SPRY domain.
Contains 1 RING-type zinc finger.

Cellular localization Cytoplasm.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors