

## Product datasheet

# Human RNF39 peptide ab45498

### Overview

---

**Product name** Human RNF39 peptide

### Description

---

**Nature** Synthetic

#### Amino Acid Sequence

**Species** Human

**Sequence** CDPRAPLRIVPAES

**Amino acids** 407 to 420

### Specifications

---

Our [Abpromise guarantee](#) covers the use of **ab45498** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**Applications** Blocking - Blocking peptide for Anti-RNF39 antibody ([ab37898](#))

**Form** Liquid

### Preparation and Storage

---

**Stability and Storage** Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.

### General Info

---

**Function** May play a role in prolonged long term-potentialiation (LTP) maintenance.

**Tissue specificity** Expressed in testis.

**Sequence similarities** Contains 1 B30.2/SPRY domain.  
Contains 1 RING-type zinc finger.

**Cellular localization** Cytoplasm.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors