

Product datasheet

Human TCF19 peptide ab23102

Overview

Product name Human TCF19 peptide

Description

Nature Synthetic

Amino Acid Sequence

Species Human

Sequence C-RSTAKAPSDTPAHE

Amino acids 346 to 359

Specifications

Our [Abpromise guarantee](#) covers the use of **ab23102** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications Blocking - Blocking peptide for Anti-TCF19 antibody ([ab5311](#))

Form Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

General Info

Function Potential trans-activating factor that could play an important role in the transcription of genes required for the later stages of cell cycle progression.

Sequence similarities Contains 1 FHA domain.
Contains 1 PHD-type zinc finger.

Developmental stage Growth regulated.

Post-translational modifications Phosphorylated upon DNA damage, probably by ATM or ATR.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors