

Product datasheet

Anti-Keratin 80 antibody ab107821

1 Image

Overview

| | |
|----------------------------|--|
| Product name | Anti-Keratin 80 antibody |
| Description | Rabbit polyclonal to Keratin 80 |
| Host species | Rabbit |
| Tested applications | Suitable for: WB |
| Species reactivity | Reacts with: Human |
| Immunogen | Synthetic peptide conjugated to KLH, corresponding to a region within C terminal amino acids 417-436 of Human Keratin 80 (NP_872313.2) |
| Positive control | NCI-H460 cell line lysate |

Properties

| | |
|-----------------------------|--|
| Form | Liquid |
| Storage instructions | Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term. |
| Storage buffer | Preservative: 0.09% Sodium Azide Constituents: PBS |
| Purity | Immunogen affinity purified |
| Purification notes | This antibody is purified through a protein A column, followed by peptide affinity purification. |
| Clonality | Polyclonal |
| Isotype | IgG |

Applications

Our [Abpromise guarantee](#) covers the use of **ab107821** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| Application | Abreviews | Notes |
|-------------|-----------|--|
| WB | | 1/100 - 1/500. Predicted molecular weight: 51 kDa. |


Target

Tissue specificity

Weakly expressed in tongue, but not skin or in any other tissues or organs examined.

Sequence similarities

Belongs to the intermediate filament family.

Images

Anti-Keratin 80 antibody (ab107821) at 1/100 dilution + NCI-H460 cell lysate at 35 µg

Predicted band size: 51 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors