

Product datasheet

Anti-MARVELD3 antibody ab118916

[2 Images](#)

Overview

Product name	Anti-MARVELD3 antibody
Description	Rabbit polyclonal to MARVELD3
Host species	Rabbit
Tested applications	Suitable for: WB, IHC-P
Species reactivity	Reacts with: Human Predicted to work with: Mouse, Rat, Rabbit, Guinea pig, Cow, Cat, Dog, Pig
Immunogen	Synthetic peptide, corresponding to a region within internal sequence amino acids 219-268 (SYFVLGAFSA SFSSGGGFGN NYSPFEGTE LEQVRQLDQQ YILRSPLIY) of Human MARVELD3 (NP_001017967.2). Run BLAST with ExPASy Run BLAST with NCBI
Positive control	WB: Jurkat cell lysate IHC-P: Human liver tissue

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Constituents: 97% PBS, 2% Sucrose
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab118916** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use a concentration of 0.25 µg/ml. Predicted molecular weight: 46 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

Application	Abreviews	Notes
-------------	-----------	-------


IHC-P Use a concentration of 4 - 8 µg/ml.

Target

Sequence similarities Contains 1 MARVEL domain.

Cellular localization Membrane.

Images




All lanes : Anti-MARVELD3 antibody (ab118916) at 0.25 µg/ml

Lane 1 : MW marker

Lane 2 : Jurkat cell lysate at 10 µg

Predicted band size: 46 kDa

Gel concentration: 12%.



ab118916, at 8 µg/ml, staining MARVELD3 in paraffin embedded Human liver tissue by Immunohistochemistry.

Cells with positive label are hepatocytes (indicated by arrows).

Magnification: 400X

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors