

## Product datasheet

# Anti-NHLRC2 antibody ab88725

[2 Images](#)

### Overview

<b>Product name</b>	Anti-NHLRC2 antibody
<b>Description</b>	Mouse polyclonal to NHLRC2
<b>Tested applications</b>	<b>Suitable for:</b> WB
<b>Species reactivity</b>	<b>Reacts with:</b> Human
<b>Immunogen</b>	Full length Human protein (AAH32598)
<b>Positive control</b>	Human pancreas lysate and NHLRC2 transfected 293T cell lysate

### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Storage buffer</b>	pH: 7.40 Constituent: PBS
<b>Purity</b>	Protein A purified
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

### Applications

Our [Abpromise guarantee](#) covers the use of **ab88725** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


Application	Abreviews	Notes
WB		1/500 - 1/1000. Predicted molecular weight: 79 kDa.

### Target

<b>Sequence similarities</b>	Contains 6 NHL repeats.
------------------------------	-------------------------

Contains 1 thioredoxin domain.


## Images



Anti-NHLRC2 antibody (ab88725) at 1/500 dilution + Human pancreas lysate at 50 µg

**Predicted band size :** 79 kDa

**Observed band size :** 45,79 kDa



**All lanes :** Anti-NHLRC2 antibody (ab88725) at 1/500 dilution

**Lane 1 :** NHLRC2 transfected 293T cell lysate

**Lane 2 :** non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

**Predicted band size :** 79 kDa

**Observed band size :** 40 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors